



# Discover InTray® DM-FungID™

Simple, accurate dermatophyte culture test — right in your clinic.

## Identify Dermatophyte Infections accurately and easily.

Complete the test in your own clinic

BioMed's InTray DM-FungID combines specimen collection, culture, and detection in one device. Combining several procedures into a single device saves time and money, while reducing exposure to collected samples.

- A single exposure system formulated to produce distinctive colony growth with a red color in the presence of growing dermatophytes.
- Easily identify typical dermatophyte characteristics both macro and microscopically.
- The medium inhibits most gram+ bacteria, gram- bacteria, yeast and saprophytic fungi.

### Convenient

- Ready-to-use diagnostic test medium: no plate preparation of medium necessary
  - Observe InTray directly on microscope stage
- Extended shelf life, easy storage: DM-FungID has a 27- month shelf life at room temperature, from the date of manufacture
- Use directly at the point-of-care: culture and identify in your own vet clinic.

### Growth Results and Microscopy

Once growth is visible, around 1 - 14 days, microscopy can be used for identification.

### Quality Control

ORGANISM	ATCC#	RESULTS
<i>T. mentagrophytes</i>	9533	Growth
<i>T. rubrum</i>	28188	Growth
<i>M. gypseum</i>	14683	Growth
<i>C. albicans</i>	60193	Inhibited
<i>A. brasiliensis</i>	16404	Inhibited
<i>S. aureus</i>	25923	Inhibited
<i>E. coli</i>	25922	Inhibited



Scan to learn more about FungID



Learn more at [biomeddiagnostics.com](http://biomeddiagnostics.com)

**WARNING:** This product can expose you to chemicals including Cycloheximide, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to [P65Warnings.ca.gov](http://P65Warnings.ca.gov).